

Artimes^{CTO}

Go Further and Faster

Semi-Compliance Coronary Balloon Catheter



**1.0mm
Balloon**

- Ideal for CTO
- 5F Compatible / 6F Guide Kissing
- Manicured and Balanced 0.016" tip
- Enhanced balloon flexibility and crossability

Artimes is a semi-compliance pre-dilatation catheter which offers numerous features for enhanced crossability, flexibility and extraordinary rewrap to facilitate the physician and procedure. The balloon sizes 1.0 and 1.25mm are designed for small and narrow lesions and for the treatment of CTOs.

Brosmed
embrace life

Balloon Specification

Balloon Diameter (mm)	1.0, 1.25, 1.5, 1.75, 2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, 4.0
Balloon Length (mm)	5, 10, 15, 20, 30
Nominal Pressure	6atm
Rated Burst Pressure	14atm
Average Burst Pressure	20atm
Compliance	Semi-compliance
Balloon Pleats	2/3
Tip Length (mm)	1.5/2.0/2.5
Distal Shaft OD	2.36F/2.55F/2.7F
Proximal Shaft OD	1.9F
Guiding Catheter Compatibility	5F Compatible/6F compatible for kissing
Working Length (cm)	140
Coating	Slider™ Hydro from tip to adjacent guide wire port (external surface); Eel™ coating for guide wire lumen

Ordering Information

Balloon Diameter	Balloon Length/Product Code				
	5mm	10mm	15mm	20mm	30mm
1.0mm	801-1005	801-1010	801-1015	801-1020	801-1030
1.25mm	801-1205	801-1210	801-1215	801-1220	801-1230
1.5mm	801-1505	801-1510	801-1515	801-1520	801-1530
1.75mm	801-1705	801-1710	801-1715	801-1720	801-1730
2.0mm	801-2005	801-2010	801-2015	801-2020	801-2030
2.25mm	801-2205	801-2210	801-2215	801-2220	801-2230
2.5mm	801-2505	801-2510	801-2515	801-2520	801-2530
2.75mm	801-2705	801-2710	801-2715	801-2720	801-2730
3.0mm	801-3005	801-3010	801-3015	801-3020	801-3030
3.25mm	801-3205	801-3210	801-3215	801-3220	801-3230
3.5mm	801-3505	801-3510	801-3515	801-3520	801-3530
4.0mm	801-4005	801-4010	801-4015	801-4020	801-4030

* Note: 1.0mm and 1.25mm balloon catheters are not available for sales in USA.



BrosMed Medical Co., Ltd

15th building, SMEs Venture Park SongShan Lake Hi-Tech Industrial Development Zone, Dongguan 523808, Guangdong, China

Office: +86 (769) 2289 2018 Fax: +86 (769) 2289 2016

E-mail: sales@brosmed.com www.brosmed.com



GRA-6001 /Rev.06