

COMPRESSAR

StrongArm™


SuperComfort™

System

Featuring the **SuperComfort™** Disc
New Innovation in Hands-Free Hemostasis™



SuperComfort Disc




"So simple and easy to use"



ADVANCED VASCULAR DYNAMICS
A Semler Technologies Company

StrongArm™ SuperComfort™ System



- The newest innovation in the CompressAR family, SuperComfort is used in the same way as the Comfort Discs.
- The SuperComfort Disc provides greater patient comfort and a new sure-grip surface.
- The V-notch enables precise positioning, easy access site cleaning, and clear visualization of the site to help achieve proper hemostasis.
- The SuperComfort Arm upgrades CompressAR Universal Stands currently in use, for quick easy conversion by CompressAR users with a simple Arm changeout.
- The robust StrongArm design provides greater stability during compression for hands-free hemostasis.
- The Release Lever works exactly the same as with the CompressAR Universal – same simple operation enables clinicians to quickly start using the CompressAR StrongArm SuperComfort System.

Continuing Excellence in Hands-Free Hemostasis™



Same reliable Comfort Disc design



Sure-grip stability

SuperComfort Disc

The StrongArm SuperComfort System builds on a tradition of excellence spanning more than 30 years. Providing greater clinician productivity and patient comfort, the CompressAR System has been used for hemostasis on millions of patients in cardiac nursing units, catheterization labs and interventional radiology labs worldwide. The StrongArm SuperComfort System provides the same cost-effective, convenient operation that clinicians have enjoyed with earlier versions of the CompressAR product line, and is used for the same wide array of hemostasis procedures.

The StrongArm SuperComfort Stand and the sterile SuperComfort Disc provide more stable femoral compression and greater patient comfort during hemostasis. Clinicians can use the SuperComfort System the same way that they've used the CompressAR Universal, to quickly start enjoying the new system's benefits. By avoiding manual compression, clinicians reduce muscular fatigue and risk of repetitive strain injury, and minimize exposure to blood during the hemostasis period.

The SuperComfort Arm upgrades CompressAR Universal Stands currently in use for quick, easy conversion by users. The StrongArm SuperComfort System is also available as a complete new CompressAR Stand. StrongArm SuperComfort Stands are available in versions which have Arms of either standard or extended length, and fit most hospital beds including the Hill-Rom TotalCare® Bed.

Studies have demonstrated early ambulation with use of external compression on diagnostic catheterization patients (references and data on file). Patient throughput in the cath lab can be increased through use of enhanced protocols and usage patterns.



ADVANCED VASCULAR DYNAMICS
A Semler Technologies Company

2326 NW Everett Street • Portland, OR 97210 USA
800.525.2555 • 503.223.2333 • FAX 503.223.8585
email info@compressar.com • www.compressar.com

COMPRESSAR
SYSTEM

CompressAR, StrongArm, StrongArm Super, SuperComfort, AVD, Advanced Vascular Dynamics, Hands-Free Hemostasis, Hill-Rom, and TotalCare are trademarks or registered trademarks of their respective owners.

M-2155-00 9/04 Copyright ©2004 Advanced Vascular Dynamics