




# BEIEN – DISPOSABLE LINEAR CUTTER

**BEIEN**



LC Series have application for Gastrointestinal Surgery, Resection Radical of Carcinoma Colorectal, Transverse Colectomy, etc. and Thoracic Surgery, including Pneumonectomy, LVRS, Lobectomy, etc.

## Features

1. Single knife in cartridge for each firing
2. Tissue positioning systems, ensure the tissues could not slide overflow
3. Balance leverage control systems, ensures parallel closure of the instrument's fork
  - Tissue positioning Point
  - Double balancing Points
  - Thickened Steel

More wider stapling lip, 1.8mm/1.9mm, ensure the consistency of tissue stapling.



## Product Information

### Linear Cutter Staplers Specification

Code	Description	Color	Tissue Type	Staple Heights
BLC60	Linear Cutter Staplers – 60mm	Blue	Medium	3.8 mm
BLC80	Linear Cutter Staplers – 80mm	Blue	Medium	3.8 mm
BLC100	Linear Cutter Staplers – 100mm	Blue	Medium	3.8 mm
GLC60	Linear Cutter Staplers – 60mm	Green	Thick	4.5 mm
GLC80	Linear Cutter Staplers – 80mm	Green	Thick	4.5 mm
GLC100	Linear Cutter Staplers – 100mm	Green	Thick	4.5 mm

### Cartridge Unit Specification

Code	Description	Color	Tissue Type	Staple Heights
BCA60	Cartridge Unit – 60mm	Blue	Medium	3.8 mm
BCA80	Cartridge Unit – 80mm	Blue	Medium	3.8 mm
BCA100	Cartridge Unit – 100mm	Blue	Medium	3.8 mm
GCA60	Cartridge Unit – 60mm	Green	Thick	4.5 mm
GCA80	Cartridge Unit – 80mm	Green	Thick	4.5 mm
GCA100	Cartridge Unit – 100mm	Green	Thick	4.5 mm