



Disposable Circular Staplers – CSS Series

CSS Series have application throughout the alimentary tract for end-to-end, end-to-side and side-to-side anastomoses.

Indications

Disposable circular staplers is intended to have applications throughout the alimentary tract for end-to-end, end-to-side and side-to-side anastomoses.

Features

1. Super-Thin anvil, more convenient for procedure.
2. Bigger inner cutting diameter, reduce anastomosis structure.
3. Rational arrangement of staples to reduce bleeding after procedure.
4. Tissue gap between anvil and cartridge, avoid tissue damage.
5. High quality titanium staples, bigger tension-resistance on anastomosis.



Product Information

Code	Description	Color	Inner DIA	Out DIA	Staple Heights
CSS17	Circular Stapler-17mm	White	10.0 mm	17.5 mm	4.8 mm
CSS21	Circular Stapler-21mm	Tan	13.0 mm	21.5 mm	4.8 mm
CSS23	Circular Stapler-23mm	Yellow	14.5 mm	23.5 mm	4.8 mm
CSS26	Circular Stapler-26mm	White	16.5 mm	26.0 mm	4.8 mm
CSS29	Circular Stapler-29mm	Blue	19.5 mm	29.0 mm	4.8 mm
CSS32	Circular Stapler-32mm	Purple	22.5 mm	32.0 mm	4.8 mm
CSS34	Circular Stapler-34mm	Green	24.5 mm	34.0 mm	4.8 mm