



Disposable Hemorrhoidal Cutter Staplers

Hemorrhoidal Cutter Staplers have application for control of rectal prolapse and hemorrhoidal disease.

Indications

Disposable Hemorrhoidal cutter staplers is intended to have applications throughout the anal canal to perform surgical treatment of hemorrhoid disease.

Features

1. Rational arrangement of staples to reduce bleeding after procedure.
2. Cutting tissue thickness adjustable according to prolapse situation.
3. Big volume of the cartridge to make sure enough space.
4. Tissue gap between anvil and cartridge, avoid tissue damage.



Product Information

Code	Description	Color	Inner DIA	Out DIA	Staple Heights
PPH32	Hemorrhoidal Cutter Stapler-32mm	Purple	23.0 mm	32.0 mm	4.0 mm
PPH34	Hemorrhoidal Cutter Stapler-34mm	Green	25.0 mm	34.0 mm	4.0 mm