



TCC Series have application for minimal invasive surgeries.

Indications

Intended to use in endoscopic surgery to set up an accessible channel during abdominal operation.

Features

- 1.Used for different size of devices, easy to operate, Maintained Pneumoperitoneum effectiveness
- 2.MRP patented magnetic design, ensure no leak and the puncture lip is clear.
- 3.Non-slip sleeve puncture lip.
- 4.No tissue for injury; Prevent puncture hernia.



Product Information

Code	Description	Inner Diameter	Puncture Length	Shaft Length
TCC-5-110	Trocar – 5mm port	5.8 mm	110 mm	193 mm
TCC-10-110	Trocar – 10mm port	11.0 mm	110 mm	197 mm
TCC-12-110	Trocar – 12mm port	13.0 mm	110 mm	199 mm
TCC-15-110	Trocar – 15mm port	16.0 mm	110 mm	203 mm